

Product datasheet for **GA103886**

Human PSMD13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMD13 CRISPRa kit - CRISPR gene activation of human proteasome 26S subunit, non-ATPase 13
Symbol:	PSMD13
UniProt ID:	Q9UNM6
Synonyms:	HSPC027; p40.5; Rpn9; S11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)