

## Product datasheet for **GA103890**

### Human PAF Receptor (PTAFR) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PTAFR CRISPRa kit - CRISPR gene activation of human platelet activating factor receptor
Symbol:	PTAFR
UniProt ID:	<a href="#">P25105</a>
Synonyms:	PAFR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)