

Product datasheet for **GA103892**

Human **TAS2R38** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R38 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 38
Symbol:	TAS2R38
UniProt ID:	P59533
Synonyms:	PTC; T2R38; T2R61
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)