


Product datasheet for **GA103895**

Human Prostaglandin D2 Receptor (PTGDR) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTGDR CRISPRa kit - CRISPR gene activation of human prostaglandin D2 receptor
Symbol:	PTGDR
Components:	GA103895G1, Prostaglandin D2 Receptor gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103895G2, Prostaglandin D2 Receptor gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103895G3, Prostaglandin D2 Receptor gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q13258
Synonyms:	AS1; ASRT1; DP; DP1; PTGDR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)