

## Product datasheet for **GA103899**

### Human **PTGER3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	PTGER3 CRISPRa kit - CRISPR gene activation of human prostaglandin E receptor 3
Symbol:	PTGER3
UniProt ID:	<a href="#">P43115</a>
Synonyms:	EP3; EP3-I; EP3-II; EP3-III; EP3-IV; EP3e; MGC27302; MGC141828; MGC141829
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)