

Product datasheet for **GA103900**

Human EP4 (PTGER4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTGER4 CRISPRa kit - CRISPR gene activation of human prostaglandin E receptor 4
Symbol:	PTGER4
UniProt ID:	P35408
Synonyms:	EP4; EP4R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)