

Product datasheet for **GA103901**

Human Prostaglandin F2 alpha Receptor (PTGFR) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTGFR CRISPRa kit - CRISPR gene activation of human prostaglandin F receptor
Symbol:	PTGFR
UniProt ID:	<u>P43088</u>
Synonyms:	FP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)