

Product datasheet for **GA103902**

Human PTGFRN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTGFRN CRISPRa kit - CRISPR gene activation of human prostaglandin F2 receptor inhibitor
Symbol:	PTGFRN
UniProt ID:	<u>Q9P2B2</u>
Synonyms:	CD9P-1; CD315; EWI-F; FPRP; SMAP-6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)