

Product datasheet for **GA103903**

Human Prostaglandin I2 Receptor (PTGIR) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTGIR CRISPRa kit - CRISPR gene activation of human prostaglandin I2 receptor
Symbol:	PTGIR
UniProt ID:	P43119
Synonyms:	IP; PRIPR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)