

Product datasheet for **GA103921**

Human PTPH1 (PTPN3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTPN3 CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase non-receptor type 3
Symbol:	PTPN3
UniProt ID:	P26045
Synonyms:	PTP-H1; PTPH1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)