

Product datasheet for **GA103922**

Human MEG1 (PTPN4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTPN4 CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase non-receptor type 4
Symbol:	PTPN4
UniProt ID:	P29074
Synonyms:	MEG; PTPMEG; PTPMEG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)