

Product datasheet for **GA103925**

Human MEG2 (PTPN9) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTPN9 CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase non-receptor type 9
Symbol:	PTPN9
UniProt ID:	P43378
Synonyms:	MEG2; PTPMEG2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)