

## Product datasheet for **GA103932**

### Human CD45 (PTPRC) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | PTPRC CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase receptor type C |
| Symbol:       | PTPRC  |
| UniProt ID:   | <a href="#">P08575</a>   |
| Synonyms:     | B220; CD45; CD45R; GP180; L-CA; LCA; LY5; T200   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                          |

#### Product images:



[View online >](#)