

Product datasheet for **GA103934**

Human PTPRCAP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTPRCAP CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase receptor type C associated protein
Symbol:	PTPRCAP
UniProt ID:	Q14761
Synonyms:	CD45-AP; LPAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)