

## Product datasheet for **GA103938**

### Human SAP1 (PTPRH) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PTPRH CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase receptor type H
Symbol:	PTPRH
UniProt ID:	<a href="#">Q9HD43</a>
Synonyms:	R-PTP-H; SAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)