

Product datasheet for **GA103942**

Human PTPRN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTPRN CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase receptor type N
Symbol:	PTPRN
UniProt ID:	Q16849
Synonyms:	IA-2; IA-2/PTP; IA2; ICA512; R-PTP-N
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)