

Product datasheet for **GA103946**

Human PTPRS activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | PTPRS CRISPRa kit - CRISPR gene activation of human protein tyrosine phosphatase receptor type S |
| Symbol: | PTPRS |
| UniProt ID: | Q13332 |
| Synonyms: | PTPSIGMA; R-PTP-S; R-PTP-sigma |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)