

## Product datasheet for **GA103955**

### Human Poliovirus Receptor (PVR) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PVR CRISPRa kit - CRISPR gene activation of human PVR cell adhesion molecule
Symbol:	PVR
UniProt ID:	<a href="#">P15151</a>
Synonyms:	CD155; HVED; Necl-5; NECL5; PVS; TAGE4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)