

## Product datasheet for **GA103974**

### Human QDPR activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	QDPR CRISPRa kit - CRISPR gene activation of human quinoid dihydropteridine reductase
Symbol:	QDPR
UniProt ID:	<u><a href="#">P09417</a></u>
Synonyms:	DHPR; PKU2; SDR33C1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)