

## Product datasheet for **GA104019**

### Human **RAP1GAP** activation kit by **CRISPRa**

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | RAP1GAP CRISPRa kit - CRISPR gene activation of human RAP1 GTPase activating protein |
| Symbol:       | RAP1GAP  |
| UniProt ID:   | <a href="#">P47736</a>   |
| Synonyms:     | RAP1GA1; RAP1GAP1; RAP1GAPII; RAPGAP   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)              |

#### Product images:



[View online >](#)