

## Product datasheet for **GA104029**

### Human Chemerin (RARRES2) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	RARRES2 CRISPRa kit - CRISPR gene activation of human retinoic acid receptor responder 2
Symbol:	RARRES2
UniProt ID:	<a href="#">Q99969</a>
Synonyms:	HP10433; TIG2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)