

Product datasheet for **GA104129**

Human ROM1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ROM1 CRISPRa kit - CRISPR gene activation of human retinal outer segment membrane protein 1
Symbol:	ROM1
Components:	GA104129G1, ROM1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104129G2, ROM1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104129G3, ROM1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q03395</u>
Synonyms:	ROM; ROSP1; RP7; TSPAN23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)