

Product datasheet for **GA104134**

Human RP9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RP9 CRISPRa kit - CRISPR gene activation of human RP9 pre-mRNA splicing factor
Symbol:	RP9
Components:	GA104134G1, RP9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104134G2, RP9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104134G3, RP9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TA86</u>
Synonyms:	PAP-1; PAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)