

Product datasheet for **GA104164**

Human MRPL23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRPL23 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein L23
Symbol:	MRPL23
Components:	GA104164G1, MRPL23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104164G2, MRPL23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104164G3, MRPL23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q16540
Synonyms:	L23MRP; RPL23; RPL23L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)