

## Product datasheet for **GA104235**

### Human CLIP170 (CLIP1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CLIP1 CRISPRa kit - CRISPR gene activation of human CAP-Gly domain containing linker protein 1
Symbol:	CLIP1
UniProt ID:	<a href="#">P30622</a>
Synonyms:	CLIP; CLIP-170; CLIP170; CYLN1; RSN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)