

Product datasheet for **GA104244**

Human Ryanodine receptor 2 (RYR2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RYR2 CRISPRa kit - CRISPR gene activation of human ryanodine receptor 2
Symbol:	RYR2
Components:	GA104244G1, Ryanodine receptor 2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA104244G2, Ryanodine receptor 2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA104244G3, Ryanodine receptor 2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q92736</u>
Synonyms:	ARVC2; ARVD2; RyR; RYR-2; VTSIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

