

Product datasheet for **GA104245**

Human Ryanodine Receptor 3 (RYR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RYR3 CRISPRa kit - CRISPR gene activation of human ryanodine receptor 3
Symbol:	RYR3
UniProt ID:	Q15413
Synonyms:	RYR-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)