

Product datasheet for **GA104249**

Human S100A1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100A1 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A1
Symbol:	S100A1
UniProt ID:	P23297
Synonyms:	S100; S100-alpha; S100A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)