

## Product datasheet for **GA104251**

### Human S100 alpha 2 (S100A2) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	S100A2 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A2
Symbol:	S100A2
UniProt ID:	<a href="#">P29034</a>
Synonyms:	CAN19; S100L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)