

## Product datasheet for **GA104253**

### Human S100A4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	S100A4 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A4
Symbol:	S100A4
UniProt ID:	<a href="#">P26447</a>
Synonyms:	18A2; 42A; CAPL; FSP1; MTS1; P9KA; PEL98
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)