

## Product datasheet for **GA104256**

### Human Psoriasin (S100A7) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	S100A7 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A7
Symbol:	S100A7
UniProt ID:	<a href="#">P31151</a>
Synonyms:	PSOR1; S100A7c
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)