

Product datasheet for **GA104257**

Human Calprotectin (S100A8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100A8 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A8
Symbol:	S100A8
UniProt ID:	<u>P05109</u>
Synonyms:	60B8AG; CAGA; CFAG; CGLA; CP-10; L1Ag; MA387; MIF; MRP8; NIF; P8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)