

Product datasheet for **GA104259**

Human **S100A10** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	S100A10 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A10
Symbol:	S100A10
UniProt ID:	P60903
Synonyms:	42C; ANX2L; ANX2LG; CAL1L; Ca[1]; CLP11; GP11; p10; P11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)