

Product datasheet for GA104260

Human S100A11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100A11 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A11
Symbol:	S100A11
UniProt ID:	P31949
Synonyms:	HEL-S-43; MLN70; S100C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)