

Product datasheet for **GA104263**

Human S100 beta (S100B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100B CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein B
Symbol:	S100B
UniProt ID:	<u>P04271</u>
Synonyms:	NEF; S100; S100-B; S100beta
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)