

Product datasheet for GA104267

Human SAA4 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | SAA4 CRISPRa kit - CRISPR gene activation of human serum amyloid A4, constitutive |
| Symbol: | SAA4 |
| UniProt ID: | P35542 |
| Synonyms: | C-SAA; CSAA |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)