

Product datasheet for **GA104270**

Human Retinal S antigen (SAG) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SAG CRISPRa kit - CRISPR gene activation of human S-antigen visual arrestin
Symbol:	SAG
UniProt ID:	P10523
Synonyms:	RP47; S-AG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)