

Product datasheet for **GA104290**

Human **SCN1A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SCN1A CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel alpha subunit 1
Symbol:	SCN1A
UniProt ID:	P35498
Synonyms:	EIEE6; FEB3; FEB3A; FHM3; GEFSP2; HBSCI; NAC1; Nav1.1; SCN1; SMEI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)