

## Product datasheet for **GA104292**

### Human **SCN2A** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	SCN2A CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel alpha subunit 2
Symbol:	SCN2A
UniProt ID:	<a href="#">Q99250</a>
Synonyms:	BFIC3; BFIS3; BFNIS; EIEE11; HBA; HBSCI; HBSCII; Na(v)1.2; NAC2; Nav1.2; SCN2A1; SCN2A2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)