

Product datasheet for **GA104294**

Human **SCN3A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SCN3A CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel alpha subunit 3
Symbol:	SCN3A
UniProt ID:	<u>Q9NY46</u>
Synonyms:	NAC3; Nav1.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)