

Product datasheet for **GA104295**

Human **SCN4A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SCN4A CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel alpha subunit 4
Symbol:	SCN4A
UniProt ID:	P35499
Synonyms:	CMS16; HOKPP2; HYKPP; HYPP; Na(V)1.4; NAC1A; Nav1.4; SkM1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)