

Product datasheet for **GA104296**

Human **SCN4B** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | SCN4B CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel beta subunit 4 |
| Symbol: | SCN4B |
| UniProt ID: | Q8IWT1 |
| Synonyms: | ATFB17; LQT10; Navbeta4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)