

Product datasheet for **GA104297**

Human Nav1.5 (SCN5A) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SCN5A CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel alpha subunit 5 |
| Symbol: | SCN5A |
| UniProt ID: | Q14524 |
| Synonyms: | CDCD2; CMD1E; CMPD2; HB1; HB2; HBBD; HH1; ICCD; IVF; LQT3; Nav1.5; PFHB1; SSS1; VF1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)