

Product datasheet for **GA104298**

Human **SCN7A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SCN7A CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel alpha subunit 7
Symbol:	SCN7A
UniProt ID:	Q01118
Synonyms:	NaG; Nav2.1; Nav2.2; SCN6A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)