

Product datasheet for **GA104301**

Human Nav1.8 (SCN10A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SCN10A CRISPRa kit - CRISPR gene activation of human sodium voltage-gated channel alpha subunit 10
Symbol:	SCN10A
UniProt ID:	Q9Y5Y9
Synonyms:	FEPS2; hPN3; Nav1.8; PN3; SNS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)