

Product datasheet for **GA104313**

Human CCL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCL3 CRISPRa kit - CRISPR gene activation of human C-C motif chemokine ligand 3
Symbol:	CCL3
UniProt ID:	<u>P10147</u>
Synonyms:	G0S19-1; LD78ALPHA; MIP-1-alpha; MIP1A; SCYA3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)