

## Product datasheet for **GA104314**

### Human LD78 beta (CCL3L1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	CCL3L1 CRISPRa kit - CRISPR gene activation of human C-C motif chemokine ligand 3 like 1
Symbol:	CCL3L1
UniProt ID:	<u><a href="#">P16619</a></u>
Synonyms:	464.2; D17S1718; G0S19-2; LD78; LD78-beta(1-70); LD78BETA; MIP1AP; SCYA3L; SCYA3L1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)