

Product datasheet for **GA104317**

Human MCP3 (CCL7) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CCL7 CRISPRa kit - CRISPR gene activation of human C-C motif chemokine ligand 7 |
| Symbol: | CCL7 |
| UniProt ID: | P80098 |
| Synonyms: | FIC; MARC; MCP-3; MCP3; NC28; SCYA6; SCYA7 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)