

Product datasheet for **GA104330**

Human Macrophage Inflammatory Protein 3 (CCL23) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CCL23 CRISPRa kit - CRISPR gene activation of human C-C motif chemokine ligand 23
Symbol:	CCL23
UniProt ID:	P55773
Synonyms:	CK-BETA-8; Ckb-8; Ckb-8-1; CKb8; hmrp-2a; MIP-3; MIP3; MPIF-1; SCYA23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)