

Product datasheet for **GA104337**

Human CX3CL1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CX3CL1 CRISPRa kit - CRISPR gene activation of human C-X3-C motif chemokine ligand 1
Symbol:	CX3CL1
UniProt ID:	P78423
Synonyms:	ABCD-3; C3Xkine; CXC3; CXC3C; fractalkine; neurotactin; NTN; NTT; SCYD1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)